### PSD '20 -- Lecture 9

**Coot**Manual refinement

#### Exercise 7

#### Coot Tutorial

- Open Coot tutorial PDF file.
   https://www2.mrc-lmb.cam.ac.uk/personal/pemsley/coot/web/tutorial/tutorial.pdf
   (link on PSD website)
- Open Coot.
- Extensions / Load tutorial model and data
- Follow along with Prof Bystroff (finish on your own)
- We will do up to the end of 4.4 Blob 1.
- Save screenshots for Blobs 1, 2, and 3. Submit as Exercise 7.


# Coot


- Extensions / tutorial model and data
- Displays
- Moving around
  - mouse
  - cursor
  - spacebar
  - goto atom
  - j
  - 0
  - mouse-left, right, middle, wheel

## Rotamers


Sidechains adopt discrete positions (rotamer), due to energetic preference.

Use Auto Fit Rotamer (refine menu) so fit rotamer to density.






# Crystallographic waters





Use Place Atom at Pointer (refine menu) to add as water


# Symmetry and crystal contacts

- Draw / cell and symmetry
- Master switch: Yes
- Symmetry by molecule: Display sphere, standard
- Show unit cell: Yes










# Exercise 7 Due Thu Nov 19

- Do Coot tutorial up to the end of Blobology
- For Blobs 1,2,3, follow the instructions. Take a screen shot at the end of each stage, showing the newly re-modeled blob.
- Paste the three images into a file and print (black/ white is fine). Add comments as desired.

#### Summary flowcharts

### From crystal to data



#### I is always a relative term


Bigger crystal, higher I

Better crystal, higher I


Longer exposure, higher I

More intense Xrays, higher I


### From data to Patterson map




### From Patterson map to phases



#### From Phases to Model




density modification?




Model

### Experimental sources of error



horizontally polarized X-rays

### Experimental sources of error




A problem for *synchrotron X-rays*. Solution: Use an external flux meter + scaling.


colimation


Wide beam means high background, large spots, spot overlap. Narrow beam means longer exposures, uneven exposure of crystal

variable wavelength


Spots may be radially smeared. Solution: Use *monochromator*.







# Computational Sources of error



Luzatti or  $\Sigma_A$  plot will estimate errors. Real-space R will locate errors.

Model data/parameter ratio

phase bias

bad geometry

Omit maps,  $2F_o$ - $F_c$  maps correct for phase bias.

PROCHECK, Molprobity find bad contacts, bad rotamers, etc. 22

## Thought experiment

What is the phase error for a 4Å resolution reflection if the coordinate error is 1Å? (in a crystal with just one atom)

Positional uncertainty in the context of the Bragg planes gives the phase uncertainty of that atom's contribution to that reflection.

If the error in atomic position is  $1\text{\AA}$ , and the Bragg plane separation is  $4\text{\AA}$ , then the error in phase is  $(1\text{\AA}/4\text{Å})*360°=90°$ 

#### Review

- systematic absences in reciprocal space are caused by lattice periodicity of the projected density in screw rotation (or glide plane)
- Sources of error are experimental and (mostly) computational
- Experimental error is reduced by scaling.
- Computational error is reduced by refinement.
- Using R-free avoids overfitting.
- R-free correlates with phase error.
- Phase error correlates with coordinate error.
- Therefore, R-free correlates with coordinate error.

### Exercise 8 -- Reading a crystallography paper

Due mon Nov 23. Answer questions. Upload as PDF.

#### Download the PDF linked to "Ex 8 paper"

**Read** the section labeled "Structure Determination"

**Explain:** "All non-hydrogen atoms were refined anisotropically."

**Read** Table 1.

**Answer**: (1) What is the meaning of the second resolution range (in parentheses)?

- (2) Why are there values in parentheses for the other items, such a R<sub>sym</sub>?
- (3) What is Wilson B?